

## **Genotyping *N. furzeri* embryos or caudal fin tissues**

### **Material**

- DirectPCR (Tail) (Viagen #102-T)
- Proteinase K solution, RNA grade (20 mg/ml, Invitrogen)

### **DNA extraction**

1. Collect samples: embryos or caudal fin clippings.
2. Prepare the digestive solution: 2ul (for fin tissue) to 5ul (for embryos) Proteinase K per 1ml DirectPCR (Tail).
3. Incubate each sample with 100 ul digestive solution at 50-55°C for 1 to 2 hours until the sample is completely dissolved. It is also okay to leave it overnight.
4. Inactivate Proteinase K by boiling at 100°C for 10 minutes (leave the cap open for the first minute so the tube does not pop during boiling)
5. Spin at high speed to precipitate debris for clear supernatant.

### **PCR and restriction enzyme verification**

1. Running 20 ul reaction is sufficient for sequencing

10 ul GoTaq Master Mix

1 ul DNA extraction

1 ul 10 uM Primer Mix (0.5ul Forward primer + 0.5ul Reverse primer)

8 ul H<sub>2</sub>O

2. Running PCR with 40 cycles to obtain more PCR product if purification step will be done by the sequencing company (usually lower yield)

v1.1: Chi-Kuo Hu, 2024/06/21